


This range is an alternative to the classic ranges called VISIT and THER used for examinations, treatments and physiotherapy. In particular, the distinguishing feature in these tables is that height adjustment occurs while the individual sections remaining completely motionless. This makes them particularly suitable for smaller environments or to make the best use of available space, as well as to maintain the position chosen for exercises once the pulley therapy system has been set up. You can opt for a simplified base frame and configure the head section with a backrest that is 68 cm long, instead of a smaller, 52 cm long section that can come with armrests. The second type of head section can have positive and negative tilt adjustment: up to $+85^\circ / -35^\circ$.

uniX Line




Pos. 3 - 4 - 5 Type of sections

111 UNIX 1



112 UNIX 2



Pos. 6 Height adjustment

D Electric, with hand control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

W Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



Model shown - code LC 111 D2 W7



Configuration example for a treatment table from the VISIT range.



For further information on options and accessories, please refer to the first pages.

TECHNICAL FEATURES

Number of sections	2
Max safety load (kg)	170
Height adjustment	Electric
Minimum/maximum table height (cm)	55-93
Mains supply	230 V 50/60 Hz 1,5A
Actuator thrust (N)	6000
IP Protection	44
Head section adjustment	Gas spring-assisted
Head section tilt - Unix 1	+70°
Head section tilt - Unix 2	+85° -35°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	70

ACCESSORIES

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET
- AC0699 PAPER ROLL HOLDER H

PADDED SUPPORTS AND ELEMENTS:

- AC0020.W? BREATHING HOLE PLUG (Unix 1 only)
- Unix 2 comes with breathing hole plug as standard
- AC0021 RPG BAR FOR COUCH (Unix 2 only)
- AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:


Refer to the first pages

upholstery color code



Pos. 3 - 4 - 5 Type of sections

121 EVO UNIX 1



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



Model shown - code LC 121 B2 W7



TECHNICAL FEATURES

Number of sections	2
Max safety load (kg) - Electric	250
Max safety load (kg) - Hydraulic	200
Height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	50-97
Mains supply	230 V 50/60 Hz 2A
Actuator thrust (N)	8000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+70°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	85

ACCESSORIES

FOR PROTECTION AND HYGIENE:
 AC0031 PAPER ROLL HOLDER F
 AC0034 STANDARD PAPER ROLL SET
 AC0035 MAXI PAPER ROLL SET
 AC0036 PROTECTIVE SHEET SET
 AC0699 PAPER ROLL HOLDER H

PADDED SUPPORTS AND ELEMENTS:
 AC0020.W? BREATHING HOLE PLUG
 AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:
 Refer to the first pages

Configuration example for a treatment table from the VISIT range.




For further information on options and accessories, please refer to the first pages.

upholstery color code




Pos. 3 - 4 - 5 Type of sections


122 EVO UNIX 2




123 EVO UNIX PLUS



124 EVO UNIX 2 TREND



125 EVO UNIX PLUS TREND



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



Model shown - code LC123B2W7



TECHNICAL FEATURES

Number of sections	2 / 4
Max safety load (kg) - Electric	250
Max safety load (kg) - Hydraulic	200
Height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	50 - 97
Mains supply	230 V 50/60 Hz 2A
Actuator thrust (N)	8000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+85° -35°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	85 (LC122; LC123) 90 (LC124; LC125)

ACCESSORIES

FOR PROTECTION AND HYGIENE:
AC0031 PAPER ROLL HOLDER F
AC0034 STANDARD PAPER ROLL SET
AC0035 MAXI PAPER ROLL SET
AC0036 PROTECTIVE SHEET SET
AC0699 PAPER ROLL HOLDER H

PADDED SUPPORTS AND ELEMENTS:
AC0021 RPG BAR FOR COUCH
AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:
 Refer to the first pages

upholstery
color code

